



Med-Rescuer® BVM Resuscitator



Med-Rescuer[®] BVM Resuscitator

- 360° Swivel Connector
- Medication port for use with MDI or syringe
- Manometer port
- Integrated PEEP connection, no adapter necessary
- Adult, Child, and Infant sizes
- Pressure release valve standard on child and infant models, optional on adult model
- All models complete with reservoir bag and O₂ tubing
- Transparent bag with positive textured grip and sensitive response
- Folds into compact carrying size
- Range of accessories including 3M™ HEPA Filter and PEEP valve
- Latex free
- Disposable

Technical Specifications

	Adult Bag	Child Bag	Infant Bag
Bag volume	1,500ml	900ml	320ml
Oxygen reservoir volume	2,900ml	2,900ml	800ml
Length	28cm	26cm	21cm
Folded length	13cm	14cm	15cm
Weight	375g	300g	200g

Connectors:	15mm I.D. and 22mm O.D.
Operation temperature:	-15°C to 50°C
Storage temperature:	-40°C to 60°C
Pressure Release Valve:	40cm H ₂ O ± 5cm H ₂ O
Product Shelf Life:	3 Years

Ordering Information

4000	Rescuer Bag Valve Mask Resuscitator, Adult size with handle
4010	Rescuer Bag Valve Mask Resuscitator, Child size with handle
4025	Rescuer Bag Valve Mask Resuscitator, Infant size with handle
4050	Rescuer Bag Valve Mask Resuscitator, Adult size with pressure release valve
4000PV	PEEP Valve
4000F	3M™ HEPA Filter attachment

Contact Information

BLS Systems Limited

1124 South Service Road West, Oakville
Ontario, Canada, L6L 5T7

Telephone: (905)-339-1069
Fax: (905)-339-1060
E-mail: sales@blssystemsltd.com
Website: <http://www.blssystemsltd.com>



Your Choice For Life