



# Med-Rescuer<sup>®</sup> BVM Resuscitator



# Med-Rescuer® BVM Resuscitator

- 360° Swivel Connector
- Medication port for use with MDI or syringe
- Manometer port
- Integrated PEEP connection, no adapter necessary
- Adult, Child, and Infant sizes
- Pressure release valve standard on child and infant models, optional on adult model
- All models complete with reservoir bag and O<sub>2</sub> tubing
- Transparent bag with positive textured grip and sensitive response
- Folds into compact carrying size
- Range of accessories including 3M™ HEPA Filter and PEEP valve
- Latex free
- Disposable

## Technical Specifications

|                         | Adult Bag | Child Bag | Infant Bag |
|-------------------------|-----------|-----------|------------|
| Bag volume              | 1,500ml   | 900ml     | 320ml      |
| Oxygen reservoir volume | 2,000ml   | 2,000ml   | 800ml      |
| Length                  | 28cm      | 26cm      | 21cm       |
| Folded length           | 13cm      | 14cm      | 15cm       |
| Weight                  | 375g      | 300g      | 200g       |

|                         |  |
|-------------------------|--|
| Connectors:             | 15mm I.D. and 22mm O.D.                      |
| Operation temperature:  | 10°C to 50°C                                 |
| Storage temperature:    | -40°C to 60°C                                |
| Pressure Release Valve: | 40cm H <sub>2</sub> O ± 5cm H <sub>2</sub> O |
| Product Shelf Life:     | 5 Years                                      |

## Ordering Information

|               |   |
|---------------|---|
| <b>4000</b>   | Rescuer Bag Valve Mask Resuscitator, Adult size with handle                 |
| <b>4010</b>   | Rescuer Bag Valve Mask Resuscitator, Child size with handle                 |
| <b>4025</b>   | Rescuer Bag Valve Mask Resuscitator, Infant size with handle                |
| <b>4050</b>   | Rescuer Bag Valve Mask Resuscitator, Adult size with pressure release valve |
| <b>4000PV</b> | PEEP Valve  |
| <b>4000F</b>  | 3M™ HEPA Filter attachment  |

## Contact Information

### BLS Systems Limited

1124 South Service Road West, Oakville  
Ontario, Canada, L6L 5T7

**Telephone:** (905)-339-1069

**Fax:** (905)-339-1060

**E-mail:** sales@blssystemsltd.com

**Website:** <http://www.blssystemsltd.com>



*Your Choice For Life*