



Rescuer[®] Pocket Resuscitator



Rescuer[®] Pocket Resuscitator

- Clear plastic, with soft inflated cuff and fold down top, designed to fit adult or child
- No assembly required
- Patented BLS One-way Valve with 3M[™] filter. Over 99% filtration efficiency (BFE/VFE)
- Anatomical shaped cushion design for air-tight fit
- Meets all clinical and technical requirements, i.e. flowrates, tidal volume and volume pressure ratios
- Not affected by extreme temperatures
- Standard 22mm fittings - can be used with any resuscitator or bag mask
- Compact carrying case - **Many colors and custom graphics available**
- 100% disposable mask and valve
- Latex free

Ordering Information

| | |
|--------------------|--|
| 2000 | Facemask with one-way valve |
| 2000-CC | Facemask with one-way valve in carrying case |
| 2000-HS | Facemask with O ₂ , headstrap and one-way valve |
| 2025-HS-CC | Facemask with O ₂ , headstrap and one-way valve in case |
| 2025-HS-CCX | Facemask with O ₂ , headstrap, one-way valve, headstrap, skin wipe and nitrile gloves in case |
| 2050 | One-way valve with filter |

Contact Information

BLS Systems Limited

1124 South Service Road West, Oakville
Ontario, Canada, L6L 5T7

Telephone: (905)-339-1069

Fax: (905)-339-1060

E-mail: sales@blssystemsltd.com

Website: <http://www.blssystemsltd.com>



Your Choice For Life