



Rescuer[®] Pocket Resuscitator



Rescuer[®] Pocket Resuscitator

- Clear plastic, with soft inflated cuff and fold down top, designed to fit adult or child
- No assembly required
- Patented BLS One-way Valve with biological filter - allows no backflow of air or contaminants from patient to rescuer - over 98% bacterial and viral filtration efficiency
- Anatomical shaped cushion design for air-tight fit
- Meets all clinical and technical requirements, i.e. flowrates, tidal volume and volume pressure ratios
- Not affected by extreme temperatures
- Standard 22mm fittings - can be used with any resuscitator or bag mask
- Compact carrying case - **Many colors and custom graphics available**
- 100% disposable mask and valve
- Latex free

Ordering Information

2000	Facemask with one-way valve
2000-CC	Facemask with one-way valve in carrying case
2000-HS	Facemask with O ₂ , headstrap and one-way valve
2025-HS-CC	Facemask with O ₂ , headstrap and one-way valve in case
2025-HS-CCX	Facemask with O ₂ , headstrap, one-way valve, headstrap, skin wipe and vinyl gloves in case
2050	One-way valve with filter

Contact Information

BLS Systems Limited

1124 South Service Road West, Oakville
Ontario, Canada, L6L 5T7

Telephone: (905)-339-1069

Fax: (905)-339-1060

E-mail: blssystemst@telus.net

Website: <http://www.blssystemsltd.com>



Your Choice For Life