

## *BLS Systems Introduces the* **FLO<sub>2</sub>MAX Filtered Oxygen Mask**



### *The Next Generation in Oxygen Therapy*


- **Deliver from 30% - 99% Oxygen; Low - Medium - High oxygen concentrations with one mask.**
- **Compatible with standard Nebulizer for Drug therapy.**
- **All patient exhalation is filtered through submicron, hydrophobic HEPA filter. Protects healthcare workers *and* patients.**

#### *Other features include*

- ❑ A closed mask breathing system with integral HEPA filter maintains respiratory isolation during oxygen and drug aerosol therapy, reducing the risk of airborne infections such as SARS and influenza.
- ❑ Over 99.9% bacterial and viral filtration efficiency.
- ❑ Extra large oxygen reservoir ensures highest possible oxygen concentration while intake valve allows patient to make-up air from atmosphere if necessary.
- ❑ 3 mask sizes (Large Adult, Small Adult and Child) ensure best fit possible.
- ❑ Cost effective; Replace Low, Med, High Oxygen concentration and nebulizer masks with one mask.
- ❑ 100% latex free and disposable.
- ❑ Patented.

# The FLO<sub>2</sub>MAX Filtered Oxygen Mask

## How it works



Product Number	Description
<b>6100</b>	Large Adult FLO <sub>2</sub> MAX Mask with O <sub>2</sub> reservoir and tubing
<b>6105</b>	Small Adult FLO <sub>2</sub> MAX Mask with O <sub>2</sub> reservoir and tubing
<b>6110</b>	Child FLO <sub>2</sub> MAX Mask with O <sub>2</sub> reservoir and tubing
<b>6250-N</b>	Rescuer Nebulizer

**BLS Systems Limited**  
 1124 South Service Road West  
 Oakville, Ontario  
 Canada L6L 5T7

Telephone: (905)-339-1069  
 Fax: (905)-339-1060  
 E-mail: [sales@blssystemsltd.com](mailto:sales@blssystemsltd.com)  
 Website: <http://www.blssystemsltd.com>