



# Rescuer<sup>®</sup> Silicone BVM Resuscitator



# Rescuer<sup>®</sup> Silicone BVM Resuscitator

- Completely autoclavable and reusable
- Textured silicone bag with positive grip and sensitive response to compliance
- Versatile 360° swivel for facemask or endotracheal tube
- Child and neonate bags with 40cm safety pressure release valve
- Clear non-rebreathing patient valve
- Folds into compact convenient size
- Full range of accessories (Adult, Child & Infant size facemasks, PEEP valve)
- Designed to meet A.S.T.M. and European standards
- Durable silicone bag and mask with polycarbonate fittings
- Complete with reservoir bag and O<sub>2</sub> tubing
- Latex free

## Technical Specifications

	Adult Bag	Child Bag	Infant Bag
Bag volume	1,500ml	900ml	320ml
Oxygen reservoir volume	2,000ml	2,000ml	800ml
Length	28cm	26cm	21cm
Folded length	13cm	14cm	15cm
Weight	375g	300g	200g
Connectors:	15mm I.D. and 22mm O.D.	Materials (All components are latex-free)	
Operation temperature:	-18°C to 50°C	Bag and mask:	Silicone
Storage temperature:	-40°C to 60°C	Valve and fittings:	Polycarbonate
		Lip and butterfly valves:	Silicone
		Oxygen reservoir and tubing:	PVC

## Ordering Information

<b>4500</b>	Silicone Reusable Bag Valve Mask Resuscitator, Adult
<b>4510</b>	Silicone Reusable Bag Valve Mask Resuscitator, Child
<b>4525</b>	Silicone Reusable Bag Valve Mask Resuscitator, Infant

## Contact Information

### BLS Systems Limited

1124 South Service Road West, Oakville  
Ontario, Canada, L6L 5T7

**Telephone:** (905)-339-1069

**Fax:** (905)-339-1060

**E-mail:** [blssystems@telus.net](mailto:blssystems@telus.net)

**Website:** <http://www.blssystemsltd.com>



*Your Choice For Life*